


# Loupe Adjustments

## Troubleshooter



Name \_\_\_\_\_

1) Describe the problem \_\_\_\_\_

2)\* Measure the focal length for each eye R: \_\_\_\_\_ L: \_\_\_\_\_

3) Measure focal length for both eyes: \_\_\_\_\_

4) The preferred working range (from the eye to the middle of the oral cavity). \_\_\_\_\_

5) Date of latest eyetest: \_\_\_\_\_

6) Is the frame adjusted correctly? \_\_\_\_\_

7) Is the angulation correct?



\* Measure the focal length **with the loupes on**; close one eye then measure from the opposite eye to the clearest point of focus. Repeat on other side.

